

# POWERED MEDICAL INSTRUMENTS

# product overview



**LARGE BONE**

**SMALL BONE**



**CARDIO THORACIC**

**NEURO/ENT/SPINAL**

**ARTHROSCOPY**

**CAST REMOVAL**

**AUTOPSY**

**deSoutter**  
MEDICAL

# ORTHOPAEDIC/TRAUMA

## LARGE BONE



*MDX-610  
with quick  
release drill  
attachment*

*KDX-600  
dedicated  
sagittal saw*

*MPX-610  
with sagittal  
saw attachment*

*TDZ-450  
TraumaDrive™  
with quick release  
drill attachment*

- Wide range of attachments with quick change locking.
- Aseptic and sterile battery, corded electric and air power options.



Ask for Large Bone  
Catalogue No. 167

Ask for TraumaDrive™  
Catalogue No. 188

# EXTREMITIES

## SMALL BONE



*ECV-221  
with bur  
attachment*

*MCZ-450  
with quick release  
drill attachment*

*EPV-221  
with sagittal saw  
attachment*

*MPZ-450  
with wire drive  
attachment*

- Wide range of attachments.
- Corded electric and air power options.
- Lever or foot control.
- Irrigation available for straight case models.



Ask for Small Bone  
Catalogue No. 170

# CARDIO THORACIC

## PRIMARY AND REVISION STERNOTOMIES



SDX-600  
Reciprocating  
saw

RDX-600  
Oscillating  
saw

DPX-170  
Oscillating  
saw

DPX-677  
Reciprocating  
saw

- Choice of battery or air power.
- Versatility of inferior or superior approach.
- Quick change sternum guard mounting.
- Twist lock reciprocating blade mounting.



Ask for Cardio Thoracic  
Catalogue No. 160

# NEURO/ENT/SPINAL

## HIGH SPEED



HST-100 straight bur head

HST-100 angle bur head

HST-100 cranio head


- 100,000 rpm rotational speed.
- Tool-free bur loading and nose piece changing.
- Extensive range of cutting tools and bur heads.
- Ultra slim and lightweight.



Ask for HST-100  
Catalogue No. 164

# ARTHROSCOPY

## SHAVER SYSTEM



Straight shaver blade

Angled shaver blade

- Automatic cutter recognition.
- One hand bi-directional shaver blade loading.
- Extra sharp, hardened blades deliver superior performance.
- Precise control via handpiece or foot control.



Ask for Arthroscopy  
Catalogue No. 193

# CAST REMOVAL

## OSCILLATING SAW SYSTEM



CSB-100  
Cordless Cast Saw

CC6  
Low voltage  
Cast Saw

CC5  
Low voltage  
Cast Saw with  
extraction

CC5A  
Cast Saw  
extraction  
System

- Low voltage, mains voltage and cordless.
- Low sound levels and reduced vibration.
- Complete range of blade types to cut all plaster materials.



Ask for Cast Saw  
Catalogue No. 113

# AUTOPSY

## OSCILLATING SAW SYSTEM



NS3  
Low voltage saw

NS170  
Air powered saw

- Low voltage or air operated.
- With or without extraction.
- Splash proof for easy cleaning.
- Variable speed for optimum cutting efficiency.



NS3A Saw and  
power module  
System



CNS3A Saw  
and extraction  
System

Ask for Autopsy Saw  
Catalogue No. 116

### De Soutter Medical Limited

Halton Brook Business Park, Weston Road, Aston Clinton  
Aylesbury, Bucks, HP22 5WF. UK  
tel: +44 (0) 1296 634 000 fax: +44 (0) 1296 634 033  
email: info@de-soutter.com web: www.de-soutter.com

### De Soutter Medical Australia Pty Ltd.

2/12-14 Apollo Drive, Hallam, Victoria 3803  
tel: 03 9702 4441 fax: 03 9702 4484  
email: australia@de-soutter.com web: www.de-soutter.com

All products listed in this brochure  
are backed by the De Soutter  
no-nonsense 12 month warranty.  
All instruments are CE marked.

**deSoutter**  
MEDICAL