

DISPOSABLES JUST GOT GREENER



ECO Pulse™

Pulse Lavage System

deSouther
MEDICAL



ECO Pulse™ green credentials

- No need to dispose of the entire handpiece, motor, batteries or tubing, just the low-cost pump kit
- Reduces clinical waste, helping to meet hospital's annual waste targets
- No battery disposal problems

60%
REDUCTION IN
CLINICAL WASTE

“
The ECO Pulse™ can be operated with any rotary handpiece, offering flexibility and cost saving ”

P20-120

Soft tissue nozzle without suction (10cm)



P20-121

Soft tissue nozzle with suction (10cm)

P20-125

Bone cleaning nozzle with suction (10cm)

P20-122

Femoral canal nozzle without suction (25 cm)


P20-123

Femoral canal nozzle with suction (25 cm)

P20-124

Femoral brush with suction (25 cm)



PUMP KITS	TYPE	ORDER CODE	SUPPLIED NOZZLES
 Soft tissue nozzle pattern	General purpose with suction	P20-700	P20-121
	Hip arthroplasty with suction	P20-701	P20-121 / P20-123
	Hip arthroplasty brush with suction	P20-702	P20-121 / P20-124
	General purpose without suction	P20-703	P20-120
	Hip arthroplasty without suction	P20-704	P20-120 / P20-122
 Bone cleaning nozzle pattern	Bone cleaning (high pressure) with suction	P20-705	P20-125
	Bone cleaning (high pressure) without suction	P20-706	P20-125
	Hip arthroplasty (high pressure) with suction	P20-707	P20-125 / P20-123
	Hip arthroplasty (high pressure) without suction	P20-708	P20-125 / P20-122

de Soutter
MEDICAL

www.de-soutter.com

De Soutter Medical USA Inc
224 Rolling Hill Road, Suite 12A
 Mooresville, NC 28117, USA
Tel: (866) 655 9040 Tel: +1(704) 655 9040 Fax: +1(704) 987 2035
Email: usa@de-soutter.com Web: www.de-soutter.com