



HIGH PERFORMANCE MOTOR OFFERS UNRIVALLED LARGE BONE REAMING AND SAWING CAPABILITY

The orthodrive* system is the latest innovation in high performance large bone power tools from De Soutter Medical. The MBQ high torque rotary handpiece provides the surgeon with ample power even for the most demanding acetabular reaming procedures. The new drill/ream handpiece produces 20% more reaming torque and the powerful dedicated sagittal saw cuts through the hardest of bone with ease and precision.

POWER

- Through continual engineering refinements, including our fourth generation hybrid motor controller and intelligent battery charging technology, we have increased power by a further 20%
- Li-lon battery technology assures class leading power to weight ratio.

VERSATILITY

- Exceptional system versatility. 15 different rotary attachments, plus 8 different function attachments, including: reciprocating and mini sagittal saws, wire/pin driver and a unique pulse lavage system.
- In-built, quickchange drill/ream gearbox eliminates dedicated reaming attachments, reducing cost and weight.
- Optional twin trigger rotary handpiece with instant forward, reverse and unique oscillating drill/tap modes.

RELIABILITY

- Proven system reliability enhanced through continuous engineering developments, including a hermetically sealed hybrid motor controller and unique intelligent battery charging technology.
- New high performance surface finish protects against the most aggressive detergents whilst offering a highly scratch resistant barrier.

Whether your preference is for battery - either aseptic or sterile - or indeed corded power, the orthodrive® system provides total OR flexibility.









MBQ MODULAR ELECTRIC SYSTEM



	Model	Order code	Drilling Speed	Reaming Speed	Maximum reaming torque	Cannulation
orthodale:	MBQ-700	1289774	0 - 1350 rpm	0 - 300 rpm	20 Nm	Ø 4.4 mm
orthodrive*	MBQ-701	1289884	0 - 1350 rpm	0 - 300 rpm	20 Nm	Ø 4.4 mm
FI H						

DB DEDICATED ELECTRIC SAW SYSTEM



orthodrive®

orthodrive

ECOPULSE TM PULSE LAVAGE SYSTEM

"The Next Generation in Pulse Lavage Sustainability"



The ECO Pulse doesn't require batteries, motors or electronics and therefore typically saves up to 55% in clinical waste, per procedure.

A smaller footprint allows less packaging to be used, this results in a reduction in valuable storage space by 2.5

PUMP KITS		
Order code	Туре	Supplied nozzle
P30-700	General purpose with suction	P30-121
P30-701	Hip arthroplasty with suction	P30-121/P30-123
P30-702	Hip arthroplasty brush with suction	P30-121/P30-124
P30-704	Hip arthroplasty without suction	P30-120/P30-122
P30-705	Bone cleaning (high pressure) with suction	P30-125
P30-707	Hip arthroplasty (high pressure) with suction	P30-123/P30-125
P30-708	Hip arthroplasty (high pressure) without suction	P30-122/P30-125
P30-709	Bone cleaning (fan spray) with suction	P30-127
P30-710	Hip arthroplasty (fan spray) with suction	P30-123/P30-127

NOZZLES			
	Order code	Туре	Useable length
	P30-125	Bone cleaning - high pressure with suction	6cm
	P30-127	Bone cleaning - fan spray with suction	6cm
	P30-120	Soft tissue without suction	10cm
< ::	P30-121	Shower spray with suction	6cm
—— () <u>(</u>	P30-122	Femoral canal without suction	25cm
	P30-123	Femoral canal with suction	25cm
The state of the s	P30-124	Femoral brush with suction	25cm

Flow rate Variable 0 - 1.0 litres/min Pulse rate Variable 0 - 1350 pulses/min * Based on AB-703 Performance figures dependent on battery type.

MBQ DRIVES AND ATTACHMENTS

	Order code	Model	Chuck description	Final drive speed	Cannulation Ø mm
DRILLING AND RE	AMING A	DAPTOR	·		
	15360	-	0.5 - 4.0 mm (5/32") keyed	As handpiece	3.1
	15370	-	0.5 - 7.4 mm (1/4") keyed	As handpiece	4.2
	15380	-	0.3 - 3.0 mm (1/8") keyless	As handpiece	3.0
	15390	-	0.3 - 7.4 mm (1/4") keyless	As handpiece	4.2
	15480	-	Small A.O. Synthes/ASIF	As handpiece	2.5
_	15490	-	Large A.O. Synthes/Protec	As handpiece	4.2
	15530	-	Hudson	As handpiece	4.2
	15540	-	Trinkle	As handpiece	4.2
	15460	-	Zimmer	As handpiece	4.2
0.0	15470	-	Hudson/Zimmer combi	As handpiece	4.2
	17130	-	Harris	As handpiece	4.2
	19330	-	DHS/DCS	As handpiece	4.2
	18190	-	Stryker Triathlon®	As handpiece	2.3
	17930	-	Aesculap Targon TX	As handpiece	4.2
SAWING ATTACHM	1ENITC				
SAWING ATTACHI	ILIVIO				
			Mini sagittal saw		
三面直路	15560	NQ-700	Uses S88 type blades	0 - 24500 cpm	-
			Reciprocating saw.		
	18280	CQ-710	Uses S84 type blades	0 - 12300 cpm	-
Elita de la companya della companya					
			Reciprocating saw.		
	18390	CQ-711	Uses S22 type blades	0 - 12300 cpm	-
RADIOLUCENT DR	ILL ADA	PTOR			
			Accepts radiolucent		
S .	15510	TQ-700	drill attachment	As handpiece	-
			driii detacriment		
-					
90° RADIOLUCEN	T DRILL.	ATTACHM	ENT		
			Uses D56 type drill bits	As handpiece	-
	Name of the last			,	
	17710				
	17710				
		_			

WIRE/PIN DRIVERS





4.0

2.0

3.2

4.0

orthodrive®

BATTERY CHARGERS



The BC-700 Battery Charger requires minimal user input, simply insert battery and the BC-700 will charge, recover and recondition it to its optimum state. **The large touch-screen display** with intuitive menu navigation, displays instant battery charge status and condition. Includes a multi-lingual display and other user settings.

The BC-700 also offers a **unique healthcheck facility**, profiling the condition of an individual battery, automatically reconditioning it, to provide optimal capacity. Alternatively, a manual healthcheck can be easily selected. The BC-700 Battery Charger accepts all De Soutter Medical batteries, ensuring total backwards compatibility.



BC-700 MULTI BAY CHARGER AND MODULE OPTIONS

Charger - Order cod	Charger module description e	Charger module Order code	For battery type	Mains power cordsets - 3 metres long.	
15700	Without charger module fitted	-	-	Please specify - Order Co	8290
16230	With 4-off CA-703 charger modules	18120	AB-703/AB-704	Euro US/Canada Australia/New Zealand Swiss	8300 8310 8320
16220	With 4-off CA-700 charger modules	15750	SB-703/SB-704		14700

The BC-700 charger modules can be quickly and easily changed to suit the batteries being charged.



SINGLE BAY CHARGERS

Model	Battery type - Order code	For battery type
BC-706	17990	SB-703/SB-704
BC-707	18000	AB-703
BC-708	18010	AB-704

LITHIUM ION ASEPTIC BATTERIES

Why choose Lithium Ion battery technology

De Soutter Medical Lithium Ion batteries offer the surgical team and sterile services a number of tangible benefits over conventional battery technology.

Key advantages are:

- Li-Ion technology improves battery 'capacity to weight' ratio when compared to NiMH or NiCad chemistries.
- Significantly extends the useable run time of the power tool.
- Eliminates 'memory effect' normally associated with other battery types.
- De Soutter Medical sterile Li-Ion technology readily withstands repeated machine washing and autoclaving.
- Unique fuel gauge incorporated within De Soutter Medical Li-Ion batteries allows the surgical team to instantly assess the battery capacity - both before and during use.



SMALL BATTERY







Aseptic housing



Aseptic shield

Model	Description	Order code	Voltage	Capacity mAh
AB-703	Battery pack	17910	10.95V	2000
AH-703	Battery housing	18080	-	-
AS-703	Battery shield	18100	-	-

LARGE BATTERY



Battery



Aseptic housing



Aseptic shield

Model	Description	Order code	Voltage	Capacity mAh
AB-704	Battery pack	17920	14,6V	2000
AH-704	Battery housing	18090	-	- //
AS-704	Battery shield	18110	-	- //



orthodrive®

LITHIUM ION STERILE BATTERIES



SMALL BATTERY



Model	Order code	Voltage	Capacity mAh
SB-703	17210	13,2V	1100

LARGE BATTERY



Model	Order code	Voltage	Capacity mAh
SB-704	17220	13,2V	2500

CORDED DIRECT POWER



STERILISATION CONTAINERS







ASEPTIC BATTERY STERILISATION SYSTEM						
Set up - optional tray contents	Standard wire container and insert tray Order code	Standard insert tray Order code				
3 off handpieces 9 off attachments 2 off aseptic battery housings and shields	18150	17960				
2 off handpieces 9 off attachments 2 off aseptic battery housings and shields	18160	17970				
STERILE BATTERY STERILISATION	TERILE BATTERY STERILISATION SYSTEM					
Set up - optional tray contents	Standard wire container and insert tray Order code	Standard insert tray Order code				
3 off handpieces						
9 off attachments2 off sterile batteries	16250	17110				
	16250	17110				

Standard wire container Order code 19010 Accepts standard insert tray

Half size wire container Order code 19000 Accepts half size insert tray

Set up - optional tray contents	Half size wire container and insert tray Order code	Half size Insert tray Order code
2 off SB-70x batteries	17390	17400

STERILE FILTERED CONTAINERS



Standard sterile filtered container Order code 19250 PTFE sterile filtered container Order code 19810 Lid gasket - Order code 670023 Accepts half size insert tray.



Half size sterile filtered container Order code 19370 Lid gasket - Order code 669903 Filter paper single use pack of 100 Order code 657243



Security tag pack of 100 Order code 657253



Sterile indicator label pack of 100 Order code 657263



PTFE filter cassette Order code 671073



UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)

De Soutter Medical Limited Halton Brook Business Park, Weston Road, Aston Clinton, Aylesbury, Bucks HP22 5WF Tel: +44 (0) 1296 634 000 Email: info@de-soutter.com

AUSTRALIA

De Soutter Medical Australia Pty Ltd 2/12-14 Apollo Drive, Hallam Victoria 3803 Tel: +61 (0) 3 9702 4441 Email: australia@de-soutter.com

AUSTRIA

De Soutter Medical Austria Hauptstrasse 122-124/Büro 01 3420 Kritzendorf Osterreich Tel: +43 (0) 2243 239000 Email: austria@de-soutter.com

BELGIUM

De Soutter Medical Belgium Bessemerstraat 14 3620 Lanaken Tel: +32 (0) 89/47 15 37 Email: belgium@de-soutter.com

FRANCE

De Soutter Medical France 1 252 Avenue du Parc des Expositions 33260 La Teste de Buch Tel: +33 (0) 5 56 54 89 36 Email: france@de-soutter.com

GERMANY

De Soutter Medical Deutschland Bahnhofstraße 4 D-66625 Nohfelden Tel: +49 (0) 68 52-99 12 46 Email: deutschland@de-soutter.com

ITALY

De Soutter Medical Italia Località Fornace SNC 27022 Casorate Primo - PV Tel: +39 (0) 2 9009 4098 Email: italy@de-soutter.com

NETHERLANDS

De Soutter Medical Nederland Gelderlandhaven 2X 3433 PG Nieuwegein Tel: +31 (0) 85 0491480 Email: nederland@de-soutter.com

UNITED STATES OF AMERICA

De Soutter Medical USA, Inc. 224 Rolling Hill Road, Suite 12A Mooresville, NC 28117 Tel: +1 (704) 655 9040 Email: Usa@de-Soutter.com



For further information on De Soutter Medical powered instruments and cast saw systems please contact the address in your area.