

deSoutter
MEDICAL



Accessories

FOR SURGICAL POWERED INSTRUMENTS

www.de-soutter.com

CONTENTS

Supplied Accessories

Chuck keys	3
Wrenches	3
Lubricating oil	3
Cleaning brush	3

Battery and Electrical Components

Aseptic batteries	4
Sterile batteries	4
Battery charger.....	5
Power supply units	5
Footswitch	5

General accessories

Twist drills - non sterile.....	6
Stem cutter.....	6
Chucks and adaptors	7

Sterilisation equipment

Sterilisation cases and insert trays.....	8
---	---




Pneumatic instrument accessories

“Wash Assist” air hoses	12
Standard air hoses	13
Wash Assist air hose parts	14
Extension hoses	15
Conversion hoses	15
Air fittings	15
Multi stage regulator	15
Irrigation pump	16
Air foot control	16
Dermatome	16
Mains supply cable	16

Sales and Service Policies

Servicing and repairs.....	17
Warranty	18





SUPPLIED ACCESSORIES

CHUCK KEYS, WRENCHES AND SPANNERS			
	Used with	Description	Order code
	Q series adaptors, DX-600/500, DU-470, DZ-450/400, DV-225/220/200	Key for 0.5 - 4 mm chuck	8780
	Q series adaptors, RX-600/500, RU-470, DX-600/500, DZ-450/400, RZ-450/400, DU-470	Key for 0.5 - 7.4 mm chuck	30062
	D32-100, MV-220	Sagittal saw blade screwdriver	2250
	D34-100	Reciprocating saw blade	2260
	AV-221	wrench	14960
	CU-470, CZ-450/400, CV-220/200, BV-220/200	Wrench	9940
	ECV-221/200, BU-470, BZ-450, BV-220	Wrench	14770
	SV-220/200, MV-220/200, GD-113	Saw blade hexagon driver	9930
	SZ-450/400, SU-470	Saw blade hexagon driver	9590
	DPX-170, NS3	12 mm/22 mm spanner	190723
WIRE GUARDS			
	D35-100, D36-100	Wire guard	10360
	For use with all wire and pin drivers	Wire guard	9570
LUBRICANTS			
	For use with Large and Small bone instruments - MPX, MPZ, DPX, DPZ, EPV, D30, HST, PLX.	25 ml bottle of mineral oil	30982
	For use with Large and Small bone instruments - MPX, MPZ, DPX, DPZ, EPV, D30, HST, PLX.	50 ml bottle of mineral oil	12570
BRUSHES			
	WZ-400/450/451, PZ-400/451, WU-470/471, PU-470, WX-500/600/610, PQ-701/702, WQ-701/707	Cleaning brush	8120
	PX-610, WX-620/630, WQ-700/710	Cleaning brush	15090

BATTERY AND ELECTRICAL COMPONENTS

For use with powered instrument handpiece types MBQ, DBK, MBU, DBC, DBR, MDX and TDZ.

ASEPTIC BATTERIES

	AB-450 Small aseptic battery 10.8V NiMH Capacity 1000 mAh Order code 14510		AH-450 Small aseptic battery housing Order code 14520		AS-450 Small aseptic battery shield Order code 14530
	AB-470 Small aseptic battery 10.8V NiMH Capacity 1000 mAh Order code 14510		AH-700 Small aseptic battery housing Order code 16140		AS-700 Small aseptic battery shield Order code 16150
	AB-600 Large aseptic battery 9.6V NiMH Capacity 2100 mAh Order code 11520		AH-610 Large aseptic battery housing Order code 14430		AS-610 Large aseptic battery shield Order code 14420
	AB-701 Medium aseptic battery 10.8V NiMH Capacity 1950 mAh Order code 15340		AH-701 Medium aseptic battery housing Order code 15670		AS-701 Medium aseptic battery shield Order code 15890
	AB-702 Large aseptic battery 14.4V NiMH Capacity 1950 mAh Order code 15350		AH-702 Large aseptic battery housing Order code 18090		AS-702 Large aseptic battery shield Order code 15900
	AB-703 Li-Ion Small aseptic battery 10.95V Capacity 2000 mAh Order code 17910		AH-703 Li-Ion Small aseptic battery housing Order code 18080		AS-703 Li-Ion Small aseptic battery shield Order code 18100
	AB-704 Li-Ion Large aseptic battery 14.6V Capacity 2000 mAh Order code 17920		AH-704 Li-Ion Large aseptic battery housing Order code 18090		AS-704 Li-Ion Large aseptic battery shield Order code 18110

STERILE BATTERIES

	SB-703 Small Li-Ion sterile battery 13.2V Capacity 1100 mAh Order code 17210		SB-704 Large Li-Ion sterile battery 13.2V Capacity 2500 mAh Order code 17220		SB-600 9.6 Vdc NiCd Capacity 1200 mAh Battery Pack Order code 11630
---	--	---	--	---	---

Battery charger, power supply, power consoles and foot switch,





BATTERY CHARGERS



BC-700 Multi bay charger

* CA-450 charging adaptor required
Order code 14540




Model	Order code	Battery Type	Charger Module	Module Order code
BC-700	15700	-	None	-
	16200	AB-450*/AB-600 SB-600	CA-702	15770 
	16210	AB-701/AB-702	CA-701	15760 
	16220	SB-703/SB-704	CA-700	15750 
	16230	AB-703/AB-704	CA-703	18120 
BC-705	17980	AB-701/AB-702	-	-
BC-706	17990	SB-703/SB-704	-	-
BC-707	18000	AB-703	-	-
BC-708	18010	AB-704	-	-



Single bay charger

POWER SUPPLY



PS-701 Corded direct power supply.
Powers handpiece types MBQ and MBU.
Power supply input 240 Vac.

Order code 17870

CM-700
Corded module
Order code 15720



PS-700 Corded direct power supply.
Powers handpiece types MBQ, MBU, MDX
and TDZ.
Power supply input 240 Vac.

Order code 15710

CM-600
Corded module
Order code 11620

POWER CONSOLES



PC-450
Used with ECV, MCZ and ACZ handpieces
Power supply input 230Vac,120Vac and 100Vac.
Order code 14970

Cordset
Order code
12850



PC-451
Used with ECV, MCZ and ACZ handpieces
Power supply with irrigation
Power supply input 230Vac,120Vac and 100Vac.
Order code 14980

FS-220
Foot switch
Order code
12830



See page 16 for 3 metre mains cable.

GENERAL ACCESSORIES

Non sterile twist drills and stem cutter


TWIST DRILLS

Diameter mm	Length mm	A.O. shank Order Code	De Soutter shank Oscillating - Order Code	A.O. equivalent Order Code
1.1	60	90252		310.11
1.5 1.5	85 110	90262 90272		310.15 310.16
2.0 2.0	100 125	90282 90292		310.19 310.21
2.5 2.5	100 195	90302	12020	310.25
2.7 2.7 2.7	135 100 195	90322 90332	12030	310.26 310.28
3.2 3.2	155 145	90342 90352		310.31 310.29
3.5 3.5	155 110	90362 90372		310.35 310.37
4.5 4.5	145 195	90382 90392		310.44 310.48
6.0	195	90402		310.60



Minimum order quantity: 5-off mixed or same type. Material: stainless steel.

STEM CUTTER



Stem cutter 38 mm dia.
Diamond disc and arbor assembly
Order code 3940

ERGOPLUS™ ACCESSORIES



Lever for ECV-220 - Order code 14870



Irrigation clip
for ECV-221
Order code
639553



Blanking ring
for ECV-221
Order code
636573

All drives on this page are cannulated excluding A.O. drill types marked *

CHUCKS AND ADAPTORS

Output drive		Input drive - Order code								
		 Q Series	 Hudson	 Zimmer	 Harris	 Jacobs	 DIN 58809	 A.O. Reamer	 Trinkle	 A.O. Drill
	Jacobs 4mm	15360								6270*
	Jacobs 7.4mm (1/4")	15370	5440	8800	8810		8820	7970	8830	11550*
	3 mm (1/8") Keyless	15380	12320	12430				12310		
	7.4 mm (1/4") Keyless	15390	12340	12440				12330		
	ISO 6360									11270*
	A.O. Drill	15480	5470	10260	10270	10280	10290	7980	10300	
	A.O. Reamer	15490	5480	10170	10180	10190	10200		10210	
	Hudson	11530		10000	10010	10020	10030	6560	10040	11580*
	Zimmer	15460	5450		10050	10060	10070	6590	10080	
	Trinkle	15540	7910	10220	10230	10240	10250	6580		
	Harris	17130	5460	10090		10100	10110	6550	10120	
	DHS	15500	11260					7770		
	DCS		17510							
	Targon TX	17930								
	DIN 58809		5490	10130	10140	10150		6570	10160	
	Biomet drill									12890*

STERILISATION EQUIPMENT


Sterilisation containers

STANDARD CONTAINER AND UNIVERSAL TRAY




Universal sterilisation case and tray order code: 9350
Universal tray order code: 631343
Lid order code: 630843
Base order code: 630833
Blue mat order code: 616153

Will accept all standard trays



WIRE CONTAINER



Wire sterilisation container and lid.
Order code 655073
Accepts both sterile and aseptic insert trays

STERILE FILTERED CONTAINER



Order code 17410
Accepts standard or wire container plus tray



Filter paper - single use
pack of 100
Order code 657243



Security tag,
pack of 100
Order code 657253

Sterile indicator label, pack of 100 Order code 657263


HALF SIZE WIRE CONTAINER



Wire Sterilisation Container and lid.
Order code 655593


Large bone battery and pneumatic sterilisation systems

ORTHODRIVE® BATTERY SYSTEM




Set up	System	Standard Container and insert tray Order code	Insert tray Order code
3 MBQ/DBK/DBC handpieces plus 9 attachments	Aseptic NiMH	17000	16990
	Aseptic Li-Ion	18150	17960
2 MBQ/DBK/DBC handpieces plus 11 attachments	Aseptic NiMH	16240	16360
	Aseptic Li-Ion	18160	17970
3 MBQ/DBK/DBC handpieces plus 9 attachments	Sterile	16250	17110

ORTHODRIVE® LITE BATTERY SYSTEM




Set up	System	Standard Container and insert tray Order code	Insert tray Order code
1 MBQ-707/708 handpiece plus 8 attachments	Aseptic NiMH	17820	17540
	Aseptic Li-Ion	18130	17900
1 MBQ-707/708 handpiece plus 8 attachments	Sterile	17830	17550

MULTIDRIVE™ BATTERY SYSTEM



Set up	System	Standard Container and insert tray Order code	Insert tray Order code
3 MDX/CDX/KDX handpieces plus 7 attachments	Aseptic	12170	12180
2 MDX/KDX handpieces plus 7 attachments	Aseptic	14440	14450
3 MDX/CDX/KDX handpieces plus 7 attachments	Sterile	11670	11690

MULTIDRIVE™ PNEUMATIC SYSTEM



Set up	Standard Container and insert tray Order code	Insert tray Order code
3 MPX/DPX handpieces	14630	14660
2 MPX/DPX handpieces	14640	14650




Washing cap for MPX, SPX and DPX - Order code 608073

Insert trays shown complete with handpieces and other items are for illustration purposes only.

STERILISATION EQUIPMENT


Trauma and Small Bone sterilisation systems

TRAUMADRIVE™ BATTERY SYSTEM




Set up	System	Standard Container and insert tray Order code	Insert tray Order code
2 MBU handpieces plus 10 attachment	Aseptic NiMH	16380	16970
	Aseptic Li ion	18140	18060
2 MBU handpieces plus 10 attachment	Sterile	16940	16980

MULTIDRIVE™ ELECTRIC - CORDED SYSTEM



Set up	Standard Container and insert tray Order code	Insert tray Order code
1 MCZ and 1 ECV handpiece, 5 ECV attachments, 8 MCZ attachments and cord set	14390	14410

MULTIDRIVE™ - PNEUMATIC SYSTEM



Set up	Standard Container and insert tray Order code	Insert tray Order code
1 MPZ and 1 EPV handpiece, 5 EPV attachments, 8 MPZ attachments and air hose	14390	14410



Washing cap for EPV, MPZ handpieces
Order code 9950

MICRODRIVE™ - PNEUMATIC SYSTEM




D30 Series Handpieces
Order code 7540



Washing cap for D30 Series and PLX handpieces
Order code 8110

STERNUDRIVE™ BATTERY SYSTEM




Model	System	Container and insert tray Order code	Insert tray only Order code
DBC-470	Aseptic and Sterile	17140	17150
DBR-700	Sterile only	17890	17880

STERNUDRIVE™ PNEUMATIC SYSTEM




Model	Container and insert tray Order code	Insert tray only Order code
DPC-470	17160	17170
DPX-170	9350*	631343

* Refer to page 8. Insert tray silicone mat - Order code 616153

 Washing cap for MPX, SPX and DPX Order code 608073

 Washing cap for DPC-470 Order code 9950


BATTERY INSERT TRAY



For SB-700 Series batteries

Order code 17400

DERMADRIVE™ ELECTRIC SYSTEM



Model	Container and insert tray Order code
GD-113 plus accessories	8570

Plastic base, plastic lid and plastic insert

NEURO/ENT PNEUMATIC SYSTEM



Model	Container and insert tray Order code
HST/DPZ	12550

Plastic base, plastic lid and plastic insert

 Washing cap for HST/DPZ handpieces. Order code 9950

 Washing cap for HST Nose only. Order code 14830

Insert trays and Containers shown complete with handpieces and other items are for illustration purposes only.









PNEUMATIC INSTRUMENT ACCESSORIES

Wash Assist air hoses

WASH ASSIST AIR HOSES



The Wash Assist co-axial hose system is recommended where machine washing of air hoses is demanded. By releasing the probe end from hose, which is retained on the hose via the rubber strap, the hose ends can be joined and sealed prior to washing.

Probe end	For model type	Complete hose length - Order code		
		3 metre	4 metre	5 metre
 MA7-Schrader	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	13000	13100	13200
		13300	13400	13500
		13600	13700	13800
 A.O.	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	13010	13110	13210
		13310	13410	13510
		13610	13710	13810
 AGA Inst.	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	13020	13120	13220
		13320	13420	13520
		13620	13720	13820
 Draeger	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	13030	13130	13230
		13330	13430	13530
		13630	13730	13830
 Heyer	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	13040	13140	13240
		13340	13440	13540
		13640	13740	13840
 Air Liquide	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	13050	13150	13250
		13350	13450	13550
		13650	13750	13850
 DIN	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	13060	13160	13260
		13360	13460	13560
		13660	13760	13860
 AGA Air	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	13080	13180	13280
		13380	13480	13580
		13680	13780	13880
Flowtech	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	13070	13170	13270
		13370	13470	13570
		13670	13770	13870
NP05	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	13090	13190	13290
		13390	13490	13590
		13690	13790	13890






To special order

Standard air hoses and wash caps

STANDARD AIR HOSES



Quick release, bayonet action hose fitting with combined coaxial air supply and exhaust hose. This system is non-preferred for machine washing.

Probe end	For model type	Complete hose length - Order code		
		3 metre	4 metre	5 metre
	MA7-Schrader MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	5300	5310	5320
		10900	10910	10920
		7500	9320	7640
	A.O. MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	5330	5340	5350
		10930	10940	10950
		8690	8700	8710
	AGA Inst. MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	5360	5370	5380
		10960	10970	10980
		8720	8730	8740
	Draeger MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	5390	5400	5410
		10990	11000	11010
		8580	8590	8600
	Heyer MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470 PLX, D30 Series	7820	7830	7840
		11020	11030	11040
		8750	8760	8770
Air Liquide	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470	10390	10400	10410
		11050	11060	11070
Flowtech	MPX, DPX, SPX, MPZ, EPV, HST, DPZ-121, DPC-470	8610	8620	8630
		11080	11090	11100

To special order

HOSE WASHING CAP



Probe end	For model type		
	MPX, DPX, SPX, Order code	MPZ, EPV, HST, DPZ-121, DPC-470 Order code	PLX, D30 Series Order code
MA7 Schrader	8950	9650	11430
A.O.	10880	9670	11450
AGA	10870	9660	11440
Draeger	10890	9680	11460
Heyer	10870	9660	11440

PNEUMATIC INSTRUMENT ACCESSORIES

Wash Assist air hose parts, adaptor without air hose, air hose without adaptor

WASH ASSIST AIR HOSE PARTS



Wash assist co-axial hose.

Hose ends connected with probe detached and retained with rubber strap.



Replacement rubber strap


Order code 632393

HOSE WITHOUT ADAPTOR



For model type	Hose length - Order code		
	3 metre	4 metre	5 metre
MPX, DPX, SPX,	632303	632313	632323
MPZ, EPV, HST, DPZ-121, DPC-470	632333	632343	632333
PLX, D30 Series	632363	632373	632383


ADAPTOR WITHOUT HOSE



Adaptor type	For model type - Order code		
	MPX	MPZ	D30
MA7	632003	632103	632203
A.O.	632013	632113	632213
AGA Inst.	632023	632123	632223
Draeger	632033	632133	632233
Heyer	632043	632143	632243
Air Liquide	632053	632153	632253
DIN	632063	632163	632263
Flowtech	632073	632173	632273
AGA Air	632083	632183	632283
NPO5	632093	632193	632293

Extension hoses, conversion hoses, convertors and air regulator.

EXTENSION HOSES




Hose length - Order code

3 metre	6 metre
5670	6120

Hose length - Order code


3 metre	5 metre
8940	10370

CONVERSION HOSES



Probe end	Female end	Order code
A.O.	MA7	5150
AGA Inst.	MA7	7510
Draeger	MA7	7520
Air Liquide	MA7	8840
1/4" NPT	MA7	6290

CONVERTORS AND AIR REGULATOR



MA7 Probe to double MA7 female
Order code 10420



MA7 Probe to A.O. female
Order code 12400

A.O. Probe to MA7 female
Order code 12410

MA7 female to 1/4" NPT female thread
Order code 5160

1/4" BSPT male to 5/8" x 27 TPI male thread
Order code 5170



Multi stage regulator with MA7 female outlet - bull nose

Inlet pressure	Outlet pressure	Air Inlet	Order code
21 - 200 bar	0 - 10 bar	5/8" BSP	5180

PNEUMATIC INSTRUMENT ACCESSORIES AND DERMATOME

Irrigation pump, air foot control, Dermatome and mains cable

IPV IRRIGATION PUMP



IPV-220
Input - 240Vac, 120Vac, 110Vac
Order code 11560

Pump sensor IPS-220
Order code 11570



MPZ Wash Assist
air hose sub assembly



Hose length - Order code

3 metre	4 metre	5 metre
632333	632343	632353

AIR FOOT PEDAL - for use with EPV-221, MPZ-451 and HST-100



Model	Order code	Foot pedal features	
		Inlet connection	Outlet connection
FV-220	14900	Integral 3 metre supply hose with MA7 probe	Accepts any MPZ "Wash Assist" air hose sub assembly - see above
FV-221	14910	Accepts Draeger Wash Assist hoses with air being exhausted at wall fitting	Accepts any MPZ "Wash Assist" air hose sub assembly - see above
FV-222	15000	Integral 3 metre supply hose with MA7 probe	MA7 outlet

DERMATOME



Model	230 V UK plug Order code	230 V Euro plug Order code	120 V US plug Order code
CM2 motor controller and footswitch	1280884	1280994	1281004



Sole plate 23 mm cutting width Order code 17340
Sole plate 48 mm cutting width Order code 17260
Sole plate 73 mm cutting width Order code 17270

MAINS CABLE



Used with:
BC-600/610/700 Battery Charger
PC-450/451 Power Console
PS-600/700/701 Power Supply
IPV-220 Irrigation Pump

UK Order code 8290
Euro Order code 8300
US/Canada Order code 8310
Australia Order code 8320
New Zealand Order code 8320
Swiss Order code 14700

Cable length 3 metres

SALES AND SERVICE POLICIES

Servicing and Repairs

Our equipment is supported by a number of strategically located servicing and repair centres. All powered instruments and attachments should be regularly serviced to ensure optimum performance and reliability. Regular servicing will minimise the risk of unexpected breakdowns. Servicing includes a comprehensive functional check, cleaning and testing, followed by the replacement of any worn or damaged components. Because of specialist techniques used in the manufacture and maintenance of our equipment, no user or third party servicing is possible.

To return equipment for repair or warranty, please ensure the following:

- Disinfect/sterilise the equipment in accordance with reprocessing instructions outlined in the products Operating Instructions.
- Record the serial number of the equipment being returned with the return documentation and enclose a brief statement describing the reason for the equipment's return together with a contact information. (contact name, email address and telephone number).
- For warranty returns, please state the original invoice number, date of purchase and a brief statement describing the reasons for the equipment's return.
- For repair/servicing, please enclose a purchase order number to cover the Minimum Repair/Inspection Charge.
- Pack instruments securely, as any in-transit damage is not covered by the warranty.
- The cost of all repairs will always be advised by De Soutter Medical and approved by the customer before any work commences.

Lithium Ion batteries - returns and transportation

Please see individual Operating Instructions for our policy for disposal, return and transportation of Lithium Ion batteries. Lithium Ion battery packs should not be returned or transported if there are any signs of damage. If in doubt, please seek guidance from your local De Soutter Medical representative.

- When returning a battery pack, our approved Lithium Ion battery label should be printed from the De Soutter Medical website and attached to the outer packaging. For more information on Lithium Ion battery transportation, please visit our support page "www.de-soutter.com/Support"

Returned goods

Requests for credit returns must be received within 60 days of the invoice date.

Equipment returned to De Soutter Medical for credit without our prior written authorisation will not be accepted. All returned goods will be subject to a re-stocking charge of 15% of the net original invoice value.

Credit cannot be given for items damaged during return shipment due to inadequate packaging. To receive a full credit, less the re-stocking charge, all products must be returned in full package quantities and any consumables must be in the original, undamaged packaging.

SALES AND SERVICE POLICIES

Warranty

De Soutter Medical, ("the Company"), warrants to the first purchaser or lessee ("Customer") that the instruments, attachments, parts and consumables manufactured by or for the Company (hereinafter collectively called "Products") are guaranteed to be free from defects in material and workmanship for the warranty period for a particular Product, beginning from the date of invoice.

New Products - Warranty Period

- Powered Surgical Handpieces and Attachments - Twelve (12) months.
- Cast and Autopsy Saws - Twelve (12) months.
- Battery Chargers and Power Supply/Control Units - Twelve (12) months.
- Pneumatic Hoses - Twelve (12) months.
- Handpiece cords and Power Cords - Twelve (12) months.
- Electric and Pneumatic Footswitches - Twelve (12) months.
- Handpiece Accessories, Bur Guards, Blade Guards and Adaptors - Twelve (12) months.
- Battery Packs - Six (6) months.
- Consumables i.e. single use Burs, Blades and Pump Kits - for such period expiring on the "use by date".

This limited warranty is restricted to repair or replacement by the Company, at its option, of any Product found to be defective during the warranty period. All warranties shall apply only to the original Customer. In no event shall the Company be liable for any anticipated profits, consequential damages or loss of time incurred by the Customer with the purchase or use of any Product.

The Company shall be under no liability in respect of any defect arising or resulting from fair wear and tear, improper or inadequate installation, use or maintenance, accidental or wilful damage, negligence, abnormal working conditions, failure to follow the Company's product operating instructions, misuse or modification or alteration or repair of the Products by unauthorised third parties or the Customer.

Repaired Products

The Company warrants to the Customer that any service or repair work performed by the Company on Products not within the warranty period shall be free from defects in workmanship and materials for a period of twelve (12) months, six (6) months for pneumatic hoses, after the date the Customer received the Product(s).

The aforementioned warranty does not affect your Statutory Rights.



UNITED STATES OF AMERICA
De Soutter Medical USA, Inc.
224 Rolling Hill Road, Suite 12A
 Mooresville, NC 28117
Tel: +1 (704) 655 9040
fax: +1 (704) 987 2035
Email: usa@de-soutter.com

UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)
De Soutter Medical Limited
Halton Brook Business Park, Weston Road
Aston Clinton, Aylesbury, Bucks HP22 5WF
Tel: +44 (0) 1296 634 000
Fax: +44 (0) 1296 634 033
Email: info@de-soutter.com

AUSTRALIA
De Soutter Medical Australia Pty Ltd
2/12-14 Apollo Drive, Hallam
Victoria 3803
Tel: +61 (0) 3 9702 4441
Fax: +61 (0) 3 9702 4484
Email: australia@de-soutter.com

AUSTRIA
De Soutter Medical Austria
Zweigniederlassung Österreich
Dietrichsteingasse 10, A-3400 Klosterneuburg
Tel: +43 (0) 676 96 71 770
Fax: +43 (0) 2243 21 656
Email: austria@de-soutter.com

BELGIUM
De Soutter Medical Belgium
In De Bruel 30
3620 Lanaken
Tel: +32 (0) 89/47 15 37
Fax: +32 (0) 89/70 12 19
Email: belgium@de-soutter.com

FRANCE
De Soutter Medical France
949 Avenue du Parc des Expositions
33260 La Teste de Buch
Tel: +33 (0) 5 56 54 89 36
Fax: +33 (0) 9 70 61 37 60
Email: france@de-soutter.com

GERMANY
De Soutter Medical Germany
Niederlassung Deutschland, Kladenfloss
D-66625 Nohfelden
Tel: +49 (0) 68 52-99 12 46
Fax: +49 (0) 68 52-99 12 47
Email: deutschland@de-soutter.com

ITALY
De Soutter Medical Italia
Località Fornace SNC
27022 Casorate Primo - PV
Tel: +39 (0) 2 9009 4098
Fax: +39 (0) 2 9009 2673
Email: italy@de-soutter.com

NETHERLANDS
De Soutter Medical Nederland
Gelderlandhaven 2X
3433 PG Nieuwegein
Tel: +31 (0) 85 0491480
Fax: +31 (0) 85 0491489
Email: nederland@de-soutter.com

deSoutter
MEDICAL

For further information on De Soutter Medical powered instruments and cast saw systems please contact the address in your area.