

# hekaDrill™




## ULTIMATE PRECISION

*The hekaDrill has been developed with innovation and precision as a guiding principle.*


Key capabilities include:

- High speed motor control - up to 85,000 rpm. Allows high torque/low speed control. (min. 200 rpm)
- Advanced motor control ensures comfort for the surgeon during extended use with minimal heat gain.
- Adjustable bur exposure lengths.
- Wireless foot pedal operates using a secure encrypted network, ensuring greater user protection.
- Adjustable acceleration and braking functionality.



*A single rotary bur allows 5 possible exposure positions offering total flexibility and reduced inventory costs.*


## ELECTRIC HANDPIECE AND ATTACHMENTS

HIGHSPEED HANDPIECE		Order code	Model	Description	Speed rpm
		103347	HDS-121	Cable length 4m	200 - 85,000
		103352	HDS-122	Lever Control	200 - 85,000
		104184	Lever		
ATTACHMENTS STRAIGHT		Order code	Model	Length - mm	Bur size - mm
		103360	S-C07	70	2.38
		103383	S-C10	100	2.38
		103394	S-C15	150	2.38
		103409	S-D09	90	3.18
		103414	S-D14	140	3.18
		103423	S-D21	210	3.18
		103438	S-D26	260	3.18
		104569	S-D-31	310	3.18
SLIMLINE					
		104035	S-C07S	70	2.38
		104056	S-C10S	100	2.38
		104312	S-C12S	120	2.38
		104062	S-C15S	150	2.38
		104429	S-C18S	180	2.38
			S-D21S	210	2.38
		104611	S-D31S	310	3.18
ANGLED					
		103446	A-C07	70	2.38
		103451	A-C10	100	2.38
		103467	A-C15	150	2.38
		103472	A-D14	140	3.18
SLIMLINE					
		104070	A-C07S	70	2.38
		104081	A-C10S	100	2.38
		104320	A-C12S	120	2.38
		104097	A-C15S	150	2.38
		104432	A-C18S	180	2.38
			A-D09S	90	3.18
			A-D21S	210	3.18
ANGLED MODULAR					
		104159	A-M01	Angled modular attachment	
FOOTED ATTACHMENT					
		103508	C-C09	90 (Small)	2.38
		103480	C-C10	100 (Medium)	2.38
		103495	C-D12	120 (Large)	3.18
METAL CUTTER		Order code	Model	Speed rpm	Chuck type
		104137	M-DO6	-	3.18
PERFORATOR ATTACHMENT					
		103519	P-P01	Hudson Used with HDS-121 at 72,000 rpm	1,000
		103686		Special order only Locking spanner for HDS-1xx supplied with handpiece	

Ask for catalogue 2677 for our range of rotary burs



Description	Model	Order code
Console - less irrigation	PC-120	103305
Console - with irrigation	PC-121	103318



The advanced hekaDrill power console is available with or without irrigation and offers numerous features:

- Simultaneous dual handpiece operation.
- Adjustable acceleration and braking functionality.
- Integrated and customisable irrigation control with anti drip mode.
- 200+ user profiles can be stored and easily accessed.
- Rings illuminate to indicate port status for ease of use.

# IRRIGATION AND FOOT PEDALS

## STERILE IRRIGATION KIT



Order code	Description	Flow Rate
S101461	C Shank	High
S101688	C Shank	Low
S101691	D Shank	-

\*Sold in Boxes of 5

## FOOT PEDALS



- Both wired and wireless foot pedal available.
- Wireless foot pedal operates using a secure encrypted network for complete protection.
- A or B tool selection.
- Foot pedal underlighting to allow for easy detection.

	Model	Order code
Wired	MFP-120	103556
Wireless	MFP-121	103562



4 metre cable for connecting and charging foot pedals via power console.

Order code: 103597

## POWER SUPPLY

Mains power cordsets -  
3 metres long.  
Please specify - Order code

UK	8290
Euro	8300
US/Canada	8310
Australia/New Zealand	8320
Swiss	14700
Brazil	21340

STERILISATION TRAYS AND CONTAINER



ELECTRIC INSTRUMENT SYSTEM

Order code: 103610  
Includes lid, container, electric tray and attachment tray.

Handpiece tray

- Optional tray contents -
- 1 off HDS-12X handpiece
  - 1 off P-P01 perforator
  - 1 off electric cable
  - 1 off spanner

Attachment tray

- Optional tray contents -
- 5 off attachments
  - 1 off spanner

CONTAINER AND TRAYS ONLY

Wire container and lid  
Order code: 103944

Tray for electric system  
Order code: 103953

Tray for attachments  
Order code: 103966



hekaDrill<sup>TM</sup>

**UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)**

De Soutter Medical Limited  
Halton Brook Business Park, Weston Road,  
Aston Clinton, Aylesbury, Bucks HP22 5WF  
Tel: +44 (0) 1296 634 000  
Email: info@de-soutter.com

**AUSTRALIA**

De Soutter Medical Australia Pty Ltd  
2/12-14 Apollo Drive, Hallam  
Victoria 3803  
Tel: +61 (0) 3 9702 4441  
Email: australia@de-soutter.com

**AUSTRIA**

De Soutter Medical Austria  
Hauptstrasse 122-124/Büro 01  
3420 Kritzensdorf  
Tel: +43 (0) 2243 239000  
Email: austria@de-soutter.com

**BELGIUM**

De Soutter Medical Belgium  
Bessemmerstraat 14  
3620 Lanaken  
Tel: +32 (0) 89/47 15 37  
Email: belgium@de-soutter.com

**FRANCE**

De Soutter Médical France  
1 252 Avenue du Parc des Expositions  
33260 La Teste de Buch  
Tél. : +33 (0) 5 56 54 89 36  
Email : france@de-soutter.com

**GERMANY**

De Soutter Medical Deutschland  
Bahnhofstraße 4  
D-66625 Nohfelden  
Tel: +49 (0) 68 52-99 12 46  
Email: deutschland@de-soutter.com

**ITALY**

De Soutter Medical Italia  
Località Fornace SNC  
27022 Casorate Primo - PV  
Tel: +39 (0) 2 9009 4098  
Email: italy@de-soutter.com

**NETHERLANDS**

De Soutter Medical Nederland  
Gelderlandhaven 2X  
3433 PG Nieuwegein  
Tel: +31 (0) 85 0491480  
Email: nederland@de-soutter.com

**UNITED STATES OF AMERICA**

De Soutter Medical USA, Inc.  
224 Rolling Hill Road, Suite 12A  
 Mooresville, NC 28117  
Tel: +1 (704) 655 9040  
Email: usa@de-soutter.com

**deSoutter**  
MEDICAL

Distributed by De Soutter Medical, manufactured by Zethon