



DPX-170S features include;

- Variable Speed Control offers maximum precision & accuracy
- High Powered pneumatic motor
- Fully Autoclavable
- Low Maintenance

### Pneumatic Saw specification

Saw model	Order code	Speed	Oscillating arc	Operating pressure
DPX-170S	224774	0 - 13000 cpm	6°	6 - 7 bar

*A range of TPLO saw blade sizes are available for purchase separately.*



### Wash Assist Hose

Model	Order code	Length	Probe
Air hose	13000	3m	MA7 Schrader
Air hose	9500	3m	1/4 NTP Connection

*Other probe connectors are available on request*



### Electric Saw specification

Saw model	Order code	Speed	Voltage
CC4S	1282774	0 - 14000 cpm	120 V
CC4S	1282884	0 - 14000 cpm	230 V

